

Expression IP5 Information Portal

ACCESSORY GUIDE

REV A










US Customers - Experience the web store at: <http://www.philips.com/newhealthcarestore>

Model 865471
June 2013



Documentation



Part Number	Description	Image
989803180291	Manual, Operator, IP5, English	Photo Not Available
989803180321	Manual, Operator, IP5, German	Photo Not Available
989803180311	Manual, Operator, IP5, French	Photo Not Available
989803180421	Manual, Operator, IP5, Spanish	Photo Not Available
989803180351	Manual, Operator, IP5, Italian	Photo Not Available
989803180281	Manual, Operator, IP5, Dutch	Photo Not Available
989803180381	Manual, Operator, IP5, Polish	Photo Not Available
989803180391	Manual, Operator, IP5, Portuguese	Photo Not Available
989803180271	Manual, Operator, IP5, Danish	Photo Not Available
989803180431	Manual, Operator, IP5, Swedish	Photo Not Available
989803180371	Manual, Operator, IP5, Norwegian	Photo Not Available
989803180401	Manual, Operator, IP5, Russian	Photo Not Available
989803180301	Manual, Operator, IP5, Finnish	Photo Not Available

Mounting Options










Part Number	Description	Image	MRI
989803181921	Desk Stand, IP5 <i>Desk Stand with tilt and height adjustment for the Expression IP5.</i>		
989803179291	Wall Arm, IP5, 16" <i>Wall arm with tilt, swivel, and height adjustment for the Expression IP5.</i>		
989803181941	Arm Extension, Strip Chart Printer <i>Extension arm to attach strip-chart printer to the 989803179291 Wall Arm for Expression IP5.</i>	Photo Not Available	
989803152951	19-inch Wall Channel <i>Allows 989803179291 Wall Arm to be mounted to a wall</i>		

Communication Accessories






Part Number	Description	Image	MRI
989803180211	Radio Kit, USB <i>Replacement 2.4GHz radio module for use with the Expression Information Portal, to communicate with wireless MRI Patient Monitors through the MRI room shielding.</i>		










989803176511	2.4GHz Flex Antenna <i>Directional antenna enhancement to improve wireless performance in certain MRI Rooms.</i>		
--------------	--	---	---

Input/output Accessories

Part Number	Description	Image	MRI
989803181881	Printer Kit, Strip Chart, USB <i>Non-MRI, USB strip-chart printer for the Expression Information Portal. Includes printer with USB cable, power adapter, and paper roll. Power cord not included.</i>		
989803152911	Printer Paper (Box of 10)		
989803178081	Bar Code Scanner, Wireless <i>MR Conditional, USB, wireless barcode scanner for the Expression Information Portal. Includes scanner and base station. Scanner can be used up to the 2,000 gauss line within the MR suite. The base station is non-MR.</i>		
989803186761	VGA Output Kit, IP5 <i>Kit to allow video output from the IP5 for viewing on an external display that can receive a VGA input. Max distance of 210ft.</i>	<p style="text-align: center;">Photo Not Available</p>	
989803178091	Keyboard and Mouse, Wireless <i>Non-MRI wireless keyboard and mouse for the Expression Information Portal. Includes keyboard, mouse, batteries, and USB wireless adapter. English language only.</i>		

Power Accessories

Part Number	Description	Image	MRI
989803179481	Power Adapter, IP5 <i>Replacement power adapter for use with the Expression IP5. Power cord not included.</i>	<p style="text-align: center;">Photo Not Available</p>	
453564177511	Universal Line Cord <i>250V cord with un terminated end</i>		
989803168211	North American Line Cord <i>Hospital grade 125V cord with US plug type</i>		

Part Number	Description	Image	MRI
453564177501	European Line Cord <i>250V cord with European (2-prong) plug type</i>		
989803173901	Brazilian Line Cord <i>250V cord with Brazilian (3-prong) plug type</i>		
989803174171	United Kingdom Line Cord <i>UK plug type</i>	Photo Not Available	
989803168221	Cord Jumper (25 feet) <i>Optional 25-foot cord jumper to extend the connection between the IP5 Power Adapter and the Line Cord</i>		
989803178101	UPS, 350VA <i>Non-MRI, 350VA uninterruptible power supply for the Expression Information Portal. Available only in US and Canada.</i>		



inspired by a better way

For more information about **Invivo Expression MRI Patient Monitoring Systems** or any of the Complete Solution products from Invivo, please visit our website at:

www.invivocorp.com or
www.ExpressionMR.com

INVIVO CORPORATION

ORLANDO, FL USA

T: 877-INVIVO1

E: INFO@INVIVOCORP.COM

WWW.INVIVOCORP.COM

WWW.EXPRESSIONMR.COM

PHILIPS

Invivo is a Philips Company

Products, services and technical specifications are subject to change without notice.

© 2013 Invivo Corporation All Rights Reserved.